

Klickmer® Assembly Protocol

Products	Klickmer® Cat. No. DX01K [fluorochrome] [size]
Recommended Use	Customization of Klickmer® with any mono-biotinylated ligand for profiling and quantitation of specific cell subsets based on receptor recognition of the ligand.
Materials Provided	Klickmer® BV421, FITC [FI], PE, APC [AP] or NONE [NO]
Materials required (not provided)	<p>Mono-biotinylated ligand to be assembled with Klickmer® must:</p> <ul style="list-style-type: none"> ▪ be prepared in aqueous buffer pH 7.0-7.5. ▪ have a biotinylation level > 75%. ▪ be free of excess biotin. <p>Dilution buffer: PBS containing 1% BSA, pH 7.</p>

Procedure Klickmer® must be assembled freshly for each experiment.

During preliminary assessment, it is recommended to test at least 3 different ratios of mono-biotinylated ligand per Klickmer® (low, medium, and high loading capacity) to optimize the reagent valency for the specific application. Too low valency may result in too low avidity to detect the target while too high valency may generate unwanted background staining.

1. Calculate the volume of assembled Klickmer® required for the experiment. *For flow cytometry experiments, see Procedural Notes for calculation example.*
2. Calculate the volume of mono-biotinylated ligand required to produce the desired volume of assembled Klickmer® using the equation below. *See Appendix 1 for calculation example.*

Your mono-biotinylated ligand volume (µL) =

$$\frac{\text{Stock Klickmer}^{\circledR} \text{ volume } (\mu\text{L}) * \text{Stock Klickmer}^{\circledR} \text{ concentration (nM)} * \text{Number of ligands per dextran}}{\text{Mono-biotinylated ligand concentration (nM)}}$$

3. In the dark, add the calculated amount of mono-biotinylated ligand into a light protected reaction tube.
4. Add Klickmer® to the mono-biotinylated ligand and mix immediately by pipetting (avoid the formation of bubbles).
5. Incubate at room temperature in the dark, for 30 min.
6. Add dilution buffer to reach the desired concentration. *Reagents needs to be pre-diluted for flow cytometry. See Procedural notes.*
7. Proceed to profiling and quantitation.

Procedural notes 1. Protocol step 6: Dilute assembled Klickmer® in aqueous buffer to a final concentration of 32 nM before staining cells for analysis using flow cytometry.

2. To analyze cells in flow cytometry using assembled Klickmer[®] reagents, follow the protocol "MHC Dextramer[®] Staining Protocol" (<https://www.immudex.com/resources/protocols/>)
3. To analyze SARS-CoV-2 Spike-specific B cells by flow cytometry using Spike Klickmer[®] reagents, follow the protocol "Application Protocol Klickmer[®] for Detection of SARS-CoV-2 Spike-specific B cells" (<https://www.immudex.com/resources/protocols/>)
4. For using Klickmer[®] reagents in flow cytometry, see the table below to see how the starting volume of Klickmer[®] correlates in test sizes:

Volume stock Klickmer [®] (160 nM)	Total volume of assembled diluted Klickmer [®] (32 nM)	Number of tests*
10 µL	50 µL	5
20 µL	100 µL	10
50 µL	250 µL	25
100 µL	500 µL	50
200 µL	1000 µL	100
1000 µL	5000 µL	500
2000 µL	10,000 µL	1000

**Each test is enough to stain 1-3 million PBMCs or 100-300 thousand clonal T cells.*

5. Always keep Klickmer[®] reagents stored at 2-8°C in the dark – the plastic vial only partially protects the reagents against light.

Technical Support

For additional Tips & Tricks, FAQs and protocols, please visit <https://www.immudex.com/resources/> or contact our support team at customer@immudex.com
Telephone: +45 3110 9292 (Denmark)

Trademarks

BV421 is equivalent to Brilliant Violet™ 421, which is a trademark or registered trademark of Becton, Dickinson and Company or its affiliates, and is used under license. Powered by BD Innovation.

Label License

Klickmer[®] BV421 is provided under an intellectual property license from Becton, Dickinson, and Company. The purchase of this product conveys to the buyer the non-transferable right to use the purchased product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components, or any materials made using the product or its components, in activity to generate revenue including (a) in manufacturing other than limited to quality control; (b) for therapeutic, diagnostic, or prophylactic purposes; or (c) for resale, regardless of whether they are resold for use in research. For information on purchasing a license to this product for purposes other than research, contact Becton, Dickinson, and Company, 10975 Torreyana Road, San Diego, California 92121 USA or bdblicensing@bd.com.

Appendix 1

Calculation Example Example to calculate the volume of mono-biotinylated ligand to assemble 20 μL Klickmer[®] (160 nM) with n mono-biotinylated ligands per dextran and a mono-biotinylated ligand concentration of (2000 nM):

$$\text{Your mono-biotinylated ligand volume } (\mu\text{L}) = \frac{(20 \mu\text{L}) * (160 \text{ nM}) * (n \text{ mono-biotinylated ligands})}{(2000 \text{ nM})}$$

Table 1. Klickmer[®] PE or BV421 Assembly (biotin acceptor sites ~20)

Ex.	Klickmer [®] vol. (μL)	Klickmer [®] conc. (nM)	n (ligands)	Ligand conc. (nM)	Ligand vol. (μL)
1	20	160	5	2000	8
2	20	160	10	2000	16
3	20	160	15	2000	24

Table 2. Klickmer[®] APC Assembly (biotin acceptor sites ~10)

Ex.	Klickmer [®] vol. (μL)	Klickmer [®] conc. (nM)	n (ligands)	Ligand conc. (nM)	Ligand vol. (μL)
1	20	160	3	2000	4.8
2	20	160	5	2000	8
3	20	160	7	2000	11.2

Table 3. Klickmer[®] FITC Assembly (biotin acceptor sites ~40)

Ex.	Klickmer [®] vol. (μL)	Klickmer [®] conc. (nM)	n (ligands)	Ligand conc. (nM)	Ligand vol. (μL)
1	20	160	10	2000	16
2	20	160	20	2000	32
3	20	160	30	2000	48