

<b>Products</b>	hCD1d/unloaded	Cat. No. XD08001M 20/100/200/1000/2000
	hCD1d/a-GalCer	Cat. No. XD08002M 20/100/200/1000/2000
	mCD1d/unloaded	Cat. No. YD08001M 20/100/200/1000/2000
	mCD1d/a-GalCer	Cat. No. YD08002M 20/100/200/1000/2000

For research use only. Not for use in diagnostic or therapeutic procedures.

**Description** The CD1d monomer complexes comprise the MR1 heavy chain,  $\beta$ 2m, loaded with or without a-GalCer. The CD1d complex is biotinylated at the C-terminus of the heavy chain. CD1d monomer is provided in liquid form in buffer containing 20 mM Tris-HCl, 50 mM NaCl, pH 8.0. The concentration of each lot will vary and be provided in the CoA.

**Sizes** 20  $\mu$ g, 100  $\mu$ g, 1 mg, 2 mg.

**Storage** Store in the dark at -80°C. Avoid repeat freeze-thaw cycles.

**Precautions** As with any product derived from biological sources, proper handling procedures should be used.

**Symbols** See [www.immudex.com/symbols](http://www.immudex.com/symbols)

**Technical support** E-mail: [customer@immudex.com](mailto:customer@immudex.com)  
Telephone: +45 3110 9292 (Denmark)

**Manufacturer** Immudex, Bredevej 2A, DK-2830 Virum, Denmark