

Products	CSS009M_x Soluble TCR Monomer is uniquely identified by its catalog number, e.g., CSS009M_x. TCR specificity is mentioned on the label e.g., HLA-A*0102 / SLLMWITQV. Affinity is mentioned on label if applicable e.g., 1 nM.
	CSS008_3 Soluble TCR Monomer is uniquely identified by its account product number, e.g., MTCRU0026. This is a custom-made product based on TCR alpha and beta chain sequences provided by the customer Soluble TCR Monomer consists of a TCR alpha and beta chain. The TCR monomer is mono-biotinylated. The product is provided in buffer containing 20 mM Tris-HCl, 50 mM NaCl, pH 8.0. For research use only. Not for use in diagnostic or therapeutic procedures
Sizes	See amount and concentration on label.
Storage	Store at -80°C. Avoid repeat freeze-thaw cycles.
Precautions	As with any product derived from biological sources, proper handling procedures should be used. For professional users.
Symbols	See www.immudex.com/symbols
Technical support	E-mail: customer@immudex.com Telephone: +45 3110 9292 (Denmark)
Manufacturer	Immudex, Bredevej 2A, DK-2830 Virum, Denmark